


Version: 1.0

**EZFAST**

# User Manual

## Blood pressure monitor

### Model: TMB-2080-K



MedM



Download MedM Health App

Thank you for selecting EZFAST Blood Pressure Monitor. Please read the user manual carefully and thoroughly so as to ensure the safe usage of this product. Keep this manual for further reference in case any issues arise.

## Table of Contents

INTRODUCTION.....	2
• General Description	
• Indications for Use	
• Contraindications	
• Measurement Principle	
• Safety Information	
• Display and Symbols	
• Name of Each Part	
• Contents/Product Includes	
• The Choice of Power Supply	
• Installing and Replacing the Batteries	
BEFORE YOU START .....	9
• Setting the Date and Time	
• Setting a user ID	
• Pair a smart device with the monitor	
• Applying the cuff	
MEASUREMENT .....	15
• Start a measurement	
• Recall the Records	
DATA MANAGEMENT.....	19
• Delete the Records	
• Tips for Measurement	
INFORMATION FOR USER .....	23
• Maintenance	
• What are systolic pressure and diastolic pressure?	
• What is the standard blood pressure classification?	
ABOUT BLOOD PRESSURE .....	25
• Irregular Heartbeat Detector	
• Why does my blood pressure fluctuate throughout the day?	
• Why do I get a different blood pressure at home compared to the hospital?	
• Is the result the same if measuring on the right arm?	
TROUBLE SHOOTING .....	28
SPECIFICATIONS.....	29
AUTHORIZED COMPONENT.....	30
EMC GUIDANCE .....	31
FCC Statement.....	34

## ♥ General Description

Thank you for selecting EZFAST arm type blood pressure monitor (TMB-2080-K). The monitor features blood pressure measurement, pulse rate measurement and the result storage.

The design provides you with two years of reliable service.

### FEATURES:

- **60 × 64 mm (2.3 × 2.5 inches) Digital LCD display**
- **Systolic, Diastolic Blood Pressure & Pulse Display**
- **Date / Time Indication**
- **Irregular Heartbeat detection**
- **2 × 250 memory storage, Guest mode available**
- **Up-to-date measuring-during-inflation technology**

## ♥ Indications for Use

The EZFAST Blood Pressure Monitor is digital monitors intended for use in measuring blood pressure and pulse rate with arm circumference ranging from 22 to 32 cm (8.6 to 12.5 inches) or 22 to 42 cm (8.6 to 16.5 inches).

It is intended for indoor, adult use only.

## ♥ Contraindications

1. The device is not suitable for use on the women who are or may be pregnant.
2. The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.














## ♥ Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a “zero pressure” equivalent to the atmospheric pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

## INTRODUCTION.

### ♥ Safety Information

The symbols below might be in the user manual, labeling or other component. They are the requirement of standard and using.

	Refer to instruction manual/ booklet To signify that the instruction manual/booklet must be read.		Type BF applied part
	Direct Current		Serial Number
	Date and Country of manufacture		Manufacturer
	Recyclable	 	The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	For indoor use only		
	Class II Equipment		
	MR Unsafe To identify an item which poses unacceptable risks to the patient, medical staff or other persons within the MR environment.		
	Caution Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.		

**CAUTION**

- \* This device is intended for indoor, home use.
- \* This device is not intended for public use.
- \* This device is portable, but it is not intended for use during patient transport.
- \* This device is not suitable for continuous monitoring during medical emergencies or operations.
- \* This device is intended for no-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm, or for any purpose other than obtaining a blood pressure measurement.
- \* This device is for adults. Do not use this device on neonates or infants. Do not use it on children unless otherwise instructed by a medical professional.
- \* Do not use on the women in pregnant, including pre-eclamptic, patients.
- \* The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.
- \* The effectiveness of this device has not been established for use:
  - on users with common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation,
  - on users with peripheral arterial disease,
  - on users undergoing intravascular therapy, or with arteriovenous (AV) shunt.
- \* Consult a medical professional before use.
- \* Do not use this device for diagnosis or treatment of any health problem or disease. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or health care professional.
- \* If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.
- \* This device may be used only for the intended use described in this manual, the manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse.
- \* Report any unexpected operation or events to the manufacturer.
- \* Do not apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.
- \* Warning: Do not kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result injury.
- \* Warning: Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- \* Warning: Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.
- \* Warning: Do not place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt, which could disrupt blood circulation and cause injuries.
- \* Do not place the cuff on the arm on the same side of a mastectomy (especially when lymph nodes have been removed). It is recommended to take measurements on the unaffected side.
- \* Do not wrap the cuff on the same arm to which another monitoring device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.
- \* Please check that the operation of the device do not result in prolonged impairment of patient blood circulation.
- \* Warning: On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.
- \* Warning: Do not use this device with high-frequency (HF) surgical equipment at the same time.


## INTRODUCTION.



### CAUTION












- \* Warning: This device is not AP/APG equipment. Do not use the device where flammable anesthetic are present, or in environments mixture with air of with oxygen or nitrous oxide.
- \* The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- \* Wireless communication equipment, such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies may cause interference that may affect the accuracy of measurements. A minimum distance of 1 foot (30 cm) should be kept from such devices during a measurement.
- \* You can use this device to take your own measurement, no third-party operator is required.
- \* Please use the device under the environment which is provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
- \* The device may require up to 30 minutes to warm up / cool down from the minimum / maximum storage temperature before it is ready for use.
- \* Warning: Excessive cuff tube lengths could cause strangulation if you don't manage them properly.
- \* Warning: Do not touch output of the batteries/adaptor and the user simultaneously.
- \* Adapter is specified as a part of ME EQUIPMENT.
- \* Warning: The power cord is considered the disconnect device for isolating this equipment from supply mains. Do not position the equipment so that it is difficult to reach or disconnect.
- \* The blood pressure monitor, its adapter, and the cuff are suitable for use within the patient environment.
- \* Warning: Do not use this device if you are allergic to polyester, nylon, or plastic.
- \* Warning: Only use accessories approved by manufacturer. Using unapproved accessories might cause damage to the unit and injure users.
- \* Warning: If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Power button immediately to release the air from the cuff.
- \* No calibration is required within two years of reliable service.
- \* Do not attempt to repair the unit yourself if it malfunctions. Only have repairs carried out by authorized service centers.
- \* At the request of authorized service personnel, circuit diagrams, component part lists, descriptions, and calibration procedures will be made available by the manufacturer or distributor.
- \* It is recommended that the performance should be checked after repair, maintenance, and every two years of use, by retesting the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50 mmHg and 200 mmHg).
- \* Warning: Do not use the device while under maintenance, or being serviced.
- \* Store your device, cuff and adapter in a clean and dry place, protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on it.
- \* Make sure the rubber tube of the cuff is not squeezed, stretched, or kinked during storage.
- \* Warning: Keep the device, cuff, and batteries away from children as they may pose a risk of choking or strangulation if used improperly.
- \* Clean both device and cuff with a soft, dry cloth. If necessary use a dampened cloth and natural detergent. Do not use alcohol, benzene, or other harsh chemicals.
- \* Do not wash the cuff in a washing machine or dishwasher!
- \* The service life of the cuff may vary by the frequency of washing, skin condition, and storage state. The typical service life is 10000 times.
- \* Dispose of accessories, detachable parts, and the device according to the local guidelines.

## ♥ Display and Symbols




- |   |                                        |    |                                                         |
|---|----------------------------------------|----|---------------------------------------------------------|
| 1 | Systolic blood pressure reading        | 8  | User ID symbol                                          |
| 2 | Blood pressure level                   | 9  | Average value symbol                                    |
| 3 | Diastolic blood pressure reading       | 10 | Data pending to transmit symbol                         |
| 4 | Pulse display                          | 11 | Bluetooth symbol                                        |
| 5 | Memory symbol                          | 12 | Cuff wearing OK symbol                                  |
| 6 | Heartbeat symbol / Irregular heartbeat | 13 | Excessive body motion detector symbol                   |
| 7 | Date / Time display                    | 14 | Battery symbol / High power symbol / Low battery symbol |

## INTRODUCTION.

SYMBOL	EXPLANATION
1	<b>Systolic blood pressure reading</b>
2	<b>Blood pressure level</b> (See the blood pressure classification on page 25-26)
3	<b>Diastolic blood pressure reading</b>
4	<b>Pulse display (Pulse rate appears after the measurement.)</b>
5	 <b>Memory symbol</b> Indicate it is in the memory mode and which the group of memory it is.
6	 <b>Heartbeat symbol</b> Flashes when detected during a measurement.
	 <b>Irregular heartbeat symbol</b> Appears when detected during a measurement. Refer to page 26 for more information.
7	Date / Time display
8	 <b>User ID symbol</b> Appears when the monitor is operated by the selected user.
9	 <b>Average value symbol</b> Appears when viewing the average value of last 3 readings under the memory mode.
10	 <b>Data pending to transmit symbol</b> Displays when the data transmission failed. The measurements are saved on the device and will be sent when a successful connection is achieved.
11	 <b>Bluetooth symbol</b> Flashes when the Bluetooth is working.
12	 <b>Cuff wearing OK symbol</b> Appears when the cuff is wrapped well.
13	 <b>Excessive body motion detector symbol</b> Appears when talking, moving, or shaking of the arm with the cuff on is detected during a measurement. NOTE: The measured blood pressure reading may not be accurate when this symbol is displayed with the reading.
14	 <b>Battery symbol / High power symbol</b> Appears when inserting the batteries or adapter. / Indicate the DC power is too high when it appears along with the error symbol «Bat H».
	 <b>Low battery symbol</b> Indicate the battery is too low when appears with the symbol «Bat L».



## ♥ Name of Each Part



## ♥ Contents/Product Includes


- Blood Pressure Monitor (TMB-2080-K)
- Cuff (Type BF applied part)  
Upper arm circumference: 22 to 42 cm (8.6 to 16.5 inches)  
Upper arm circumference: 22 to 32 cm (8.6 to 12.5 inches) (optional)
- User manual
- 4x AAA batteries (optional)
- USB cable (optional)
- AC Adapter (optional)  
Reference model: BLJ06L050100U-U

## BEFORE YOU START

### ♥ The Choice of Power Supply

1. Battery powered mode:  
6V DC 4× AAA batteries
2. AC adapter powered mode: 5V  $\equiv$  1A  
(Please only use AC adapter authorized by the manufacturer!)

Please unplug the adapter to depart from the using utility power, when you finish the measurement.




#### CAUTION

In order to get the best effect and protect your monitor, please use the right batteries and special power adapter which complies with local safety standard.



### ♥ Installing and Replacing the Batteries

- Slide off the battery cover.
- Install or replace 4 AAA size batteries according to the polarity indications inside the battery compartment.
- Place back the battery cover.

Any time the battery is low, it will display the icon  "bAt Lo" and then power off automatically after about 5 seconds.

Replace the batteries:

- When the low battery symbol appears.(see the LCD display on the right picture)
- When any button is pressed and nothing is displayed on the screen.



Low Battery






#### CAUTION


- Do not use new and used batteries together.
- Do not use different types of batteries together.
- Do not dispose the batteries in fire. Batteries may explode or leak.
- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old batteries from the device following your local recycling guidelines.

## ♥ Setting the Date and Time


It is important to set the date and time before using your blood pressure monitor for the first time, so that a correct time stamp can be assigned to each record that is stored in the memory.

(The setting range of the year: 2021—2050, Time format: 24H / 12H)







1. When the monitor is off, press and hold the  button to display the date format. Press the  or  button to switch the date format between [month/day/year] and [day/month/year].




[month/day/year] format




[day/month/year] format

2. Press the  button to confirm the date format, then the year will flash. Press the  or  button to change the year. Press and hold the  button to quickly advance the years. Press and hold the  button to quickly go backwards through the years.
3. Press the  button to confirm the year, then the month will flash. Repeat the same steps to change the month, day, time format, hour, and minute.







## BEFORE YOU START



4. After confirming the minute, the LCD will display “do nE” and the monitor will shut off after several seconds.





## ♥ Setting a user ID


There are 3 user ID    available. The  user and , each with 250 memory spaces, are designed for 2 different people to save the measured values separately. The user , no memory space, is reserved for guest.

1. When the monitor is off, press the  button to display the current user ID and blink. Press the  button to switch the user ID between the user  and .

### Note

Only user  and  can be selected in this setting mode.



2. Press the  button to confirm the selected user ID, the User ID will not blink any more and the monitor will turn off after several seconds.



## BEFORE YOU START


### ♥ Pair a smart device with the monitor


You are the intended operator of this blood pressure monitor. You can measure your blood pressure and then save and send measurement data to a smart device (such as smartphone or tablet) with Bluetooth wireless connectivity.



Download MedM Health App

- 1. Turn on Bluetooth and the app on your smart device.**  
Make sure both are ON when pair-up is proceeding.
- 2. When the monitor is off, press and hold the  button to start pair-up, the bluetooth symbol  will flash.**



- 3. If successful, the bluetooth symbol  will not flash any more and the monitor will display "do nE".**

Press the  button to turn off the device, or the monitor will automatically shut off after about 1 minute.

#### **Note**

The date and time on your monitor will automatically be set after paired with your smart device successfully.

If **unsuccessful** within 60 seconds, it is judged timeout and the monitor will shut off.

**Bluetooth Module information:**

Bluetooth Module No.: B2084

RF Frequency Range: 2402 MHz to 2480 MHz

Output Power Range:  $\leq -2.5$  dBm

Supply Voltage: 1.8-3.6 V

Transmitting Distance: 10 meters

**List of compatible devices:**

For iOS devices:

The operating system must be iOS 13.0 or more.

For Android devices:

The operating system must be Android 5.0 or more.

**CAUTION**

- Interference may occur in the vicinity of equipment marked with the following symbol . And TMB-2080-K may interfere the vicinity electrical equipment.
- Sensitive people, including pregnant women, including pre-eclamptic, patients and those patients with implanted, electrical devices should avoid using this product whenever possible.
- Keep the monitor at least 20 cm (7.8 inches) away from the human body (especially the head) when the data transmission is proceeding after the measurement.
- To enable the data transmission function, this product should be paired to Bluetooth end at 2.4 GHz.

**How to mitigate possible interference?**

- 1.The range between the device and BT end should be reasonably close, from 1 to 10 m (3.2 to 32.8 ft). Please ensure no obstacles between the device and BT end so as to obtain quality connection and to lower the RF output range.
- 2.To avoid interference, other electronic devices (particularly those with wireless transmission / Transmitter) should be kept at least 1 meter away from the monitor. And BT end so as to obtain quality connection and to lower the RF output range.

## MEASUREMENT

### ♥ Applying the cuff

1. Remove all jewelry, such as watches and bracelets from your left arm.

**Note:** If your doctor has diagnosed you with poor circulation in your left arm, use your right arm.

2. Roll or push up your sleeve to expose the skin. Make sure your sleeve is not too tight.

3. Hold your arm with your palm facing up and tie the cuff on your upper arm, align the Artery indicator  $\phi$  with the main Artery (on the inside of your arm).

**Note:** Locate the main Artery by pressing with 2 fingers approximately 2 cm (0.7 inches) above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest, that is your main Artery!

4. Make sure the bottom edge of the arm cuff is 2 to 3 cm (0.7 to 1.1 inches) above the inside elbow. Then wrap the cuff securely.




**Note:** The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

5. Sit upright in a comfortable chair with your back against the backrest of the chair. Keep your feet flat and your legs uncrossed. Place your arm resting comfortably on a flat table. The cuff worn on your arm should be placed at the same level as your right atrium of the heart.

6. Take 5-6 deep breaths and let's start measuring!

#### Helpful tips:







- Take the measurement in a silent room.
- Rest for 5 minutes before a measurement.
- Wait at least 3 minutes before another measurement. This allows your blood circulation to recover.
- Be relax and do not move and talk during the measurement procedure.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.






## ♥ Start a measurement


You can use your monitor without pairing to a smart device. To pair your monitor with a smart device refer to the foregoing chapters.

1. When the monitor is off, press the  button, the current user ID will flash.
2. Press the  button to switch the user  or  or . Press  the button to confirm the selected user ID.
3. Remain still and do not talk until the full measurement is complete. (Example shown below for User 1)

Select user ID




Cuff wrap detection



Inflating and measuring  
(a heartbeat is detected)







Display and save  
the measurement result




## MEASUREMENT



### Note

- When the selected user has history memory record, the monitor will display its last measurement record first before starting “cuff wrap detection”.
- Any time, to stop the measurement, press the  button.



4. Press the  button to turn off the monitor, or it will shut off after about 1 minute.
5. If your monitor is already paired with your smart device and both Bluetooth and app are ON, when the measurement completed, the measurement result will start transmitting. (only user 1 and 2 available) If successful, both symbols   will disappear and the monitor will display “do nE”.

Press the  button to turn off the device, or the monitor will automatically shut off after about 1 minute.






If unsuccessful within 60 seconds, it is judged timeout and the monitor will shut off. In the case of a data transmission failure, up to 250 measurements are saved on the device for each user  and  and will be sent when a successful connection is achieved.

### Note


- Both user  and  can store maximum 250 groups of record. When you pass that limit, every time you take the measurement, the monitor will prompt “FULL” first and the oldest record will be dropped from the list after the measurement.

## 6. About the irregular heartbeat and excessive body motion during the measurement.








- During a measurement, If an irregular heartbeat is detected, the  symbol will display in the measurement result. See page 26 for more information.
- During a measurement, when the excessive body motion, especially of the arm the cuff is worn is detected, the symbol  will flash about 5 seconds and detect again. If it is no longer detected, the symbol will disappear; If still detected, the symbol  will final display in the measurement result.


### Note

The measured blood pressure reading may not be accurate if this symbol is displayed.






### ♥ Recall the Records

1. When the monitor is off, press the  button, the LCD will display the total groups of the memory records, with the user ID blinking.
2. Press the  button to switch the user ID between user  and . press the  button to confirm the selected user ID.
3. Then the LCD will display the average value of last 3 readings. When the records are less than three groups, it will display the latest record. (Example shown below for User 1)
4. Press the  or  button to display the next record.





#### Note


- When your monitor is paired with your smart device and both Bluetooth and app are ON, the records will start transmitting with the symbol  blinking. Once successful, both symbols   will disappear.
- The latest record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (250) will be dropped from the list.

## ♥ Delete the Records




Delete the records by following the steps below.

**A: Delete one record (Example shown below for User 1)**

1. Enter the memory recall mode as described in the previous section «Recall the records», find out the record you want to delete.
2. Press and hold the  button about 3 seconds, the LCD will display “dEL yES” and blink.
3. Press the  button to confirm the deletion, the LCD will display “dEL donE” and then the previous record will be displayed.







### Note

Before you confirm the deletion, you could press the  or  button to switch the LCD display from “dEL yES” to “dEL no”. Then press the  button, you could stop clearing the memory.




## DATA MANAGEMENT

### B: Delete all records (Example shown below for User 1)


1. Enter the memory recall mode as described in the previous section «Recall the records».
2. Press and hold the  and  button about 5 seconds, the LCD will display “dEL ALL” and blink.
3. Press the  button to confirm the deletion, the LCD will display “dEL ALL donE”.




#### Note

Before you confirm the deletion, you could press the  or  button to switch the LCD display from “dEL ALL” to “dEL no”. Then press the  button, you could stop clearing the memory.


4. Once deleted, your readings cannot be restored. The LCD will display "--" like the following picture.

Press the  button to turn off the monitor, otherwise it will power off automatically after about 1 minute.



## ♥ Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances.



Within 1 hour after  
dinner or drinking



Immediate measurement  
after tea, coffee, smoking



Within 20 minutes after taking a bath



When talking or moving your fingers



In a very cold environment



When you want to discharge urine





## ♥ Maintenance

In order to get the best performance, please follow the instructions below.



Put in a dry place and avoid the sunshine



Avoid touching water, clean it with a dry cloth in case.



Avoid intense shaking and collisions



Avoid dusty and unstable temperature environment



Using wet cloths to remove dirt





Do not attempt to clean the reusable cuff with water and never immerse the cuff in water.

## ABOUT BLOOD PRESSURE

### ♥ What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.



### ♥ What is the standard blood pressure classification?

The chart on the right is the standard blood pressure classification published by American Heart Association (AHA).

Blood Pressure Category	Systolic mmHg (upper#)		Diastolic mmHg (lower#)
Normal	less than 120	and	less than 80
Elevated	120-129	and	less than 80
High Blood Pressure (Hypertension) Stage 1	130-139	or	80-89
High Blood Pressure (Hypertension) Stage 2	140 or higher	or	90 or higher
Hypertensive Crisis (Consult your doctor immediately)	Higher than 180	and /or	Higher than 120

**CAUTION**

Please consult a physician if your measuring result falls outside the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

## ♥ Irregular Heartbeat Detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, the monitor records all the pulse intervals and calculate the average; if there are two or more pulse intervals, the difference between each interval and the average is more than the average value of  $\pm 25\%$ , or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of  $\pm 15\%$ , the irregular heartbeat symbol appears on the display when the measurement results are appeared.

**CAUTION**

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heart-beat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

### ♥ Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times everyday. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
2. If the person takes medicine, the pressure will vary more.
3. Wait at least 3 minutes for another measurement.



### ♥ Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc. Also, there is the “white coat” effect, which means blood pressure usually increases in clinical settings.

What you need to pay attention to when you measure your blood pressure at home:

If the cuff is tied properly.

If the cuff is too tight or too loose.

If the cuff is tied on the upper arm.


If you feel anxious.

Taking 2-3 deep breaths before beginning will be better for measuring.

Advice: Relax for 4-5 minutes until you calm down.



### ♥ Is the result the same if measuring on the right arm?

It is ok for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.



## TROUBLE SHOOTING

If any abnormality arises during use, please check the following points:

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
<b>No power</b>	Display can not light up.	Batteries are depleted.	Replace with new batteries.
		Batteries are inserted incorrectly.	Insert the batteries correctly.
		Adapter is inserted incorrectly.	Insert the AC adapter correctly.
<b>High Battery</b>	bAt H shows	The battery is too high.	Replace with new batteries.
<b>Low Battery</b>	bAt Lo &  shows	The battery is too low.	Replace with new batteries.
<b>Error message</b>	E 1 shows	The cuff is not wrapped or wrapped incorrectly, or the cuff air plug is loose.	Refasten the cuff and insert air tube plug correctly then measure again.
	E 2 or  shows	Excessive body motion (such as shaking of the arm with the cuff on) or weak Pulse is detected.	Relax for 5 minutes. and then keep still, measure again.
	E 3 shows	Pulse is not detected during measuring.	Loosen the clothing on the arm and measure again.
	E 4 shows	The measurement failed.	Relax for 5 minutes and measure again.
	EEx shows	Hardware error (X can be some digital symbol, such as 1, 2, 3, etc.)	Turn off monitor and measure again. If EEx still appears on the display, please contact the retailer or our customer service.
	Err & Usb shows	Adapter error	Replace with the authorized adapter.
<b>Warning message</b>	out shows	Out of measurement range	Relax for a moment and then measure again. If the problem persists, contact your physician.

**NOTE:** If the product still does not work, contact Customer Service. Under no circumstance should you disassemble or attempt to repair the unit by yourself.


## SPECIFICATIONS

<b>External dimensions</b>	Approx.100 mm × 129 mm × 45 mm (3.9 x 5.0 x 1.7 inches)
<b>Display mode</b>	Digital LCD V.A.60 mm × 64 mm (2.3 x 2.5 inches)
<b>Weight</b>	Approx.216 g (0.4 lb) (Excluding the batteries and cuff)
<b>Measurement mode</b>	Oscillographic testing mode
<b>Mode of operation</b>	Continuous operation
<b>Measurement range</b>	Rated cuff pressure: 0 mmHg~299 mmHg Measurement pressure: SYS: 60 mmHg ~ 230 mmHg DIA: 40 mmHg ~ 130 mmHg Pulse value: (40-199) beat/minute
<b>Accuracy</b>	Pressure: 5 °C-40 °C within ±3 mmHg Pulse value: ±5%
<b>Normal working condition</b>	A temperature range of: 5 to 40 °C (41 to 104 °F) A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa An atmospheric pressure range of : 700 hPa to 1060 hPa
<b>Storage &amp; transportation condition</b>	Temperature: -20 to 60 °C (-4 to 140 °F) A relative humidity range of ≤ 93%, non-condensing, at a water vapour pressure up to 50 hPa
<b>Measurement perimeter of the upper arm</b>	About 22 to 32 cm (8.6 to 12.5 inches) or 22 to 42 cm (8.6 to 16.5 inches)
<b>Degree of protection</b>	Type BF applied part
<b>Protection against ingress of water</b>	IP21 It means the device could protected against solid foreign objects of 12.5 mm (0.4 inches) and greater, and protect against vertically falling water drops.
<b>Device Classification</b>	Battery Powered Mode: Internally Powered ME Equipment AC Adapter Powered Mode: Class II ME Equipment

**WARNING:** No modification of this equipment is allowed.

## ♥ Authorized Component

please use the authorized adapter (optional).



### Adapter

Type: BLJ06L050100U-U

Input: 100-240V, 50-60Hz, 0.2A max

Output: 5V  $\equiv$  1000 mA

## ♥ Contact Information

### i-SENS USA, Inc.

3545 Lomita Blvd., Unit D, Torrance, CA 90505 USA

Customer Service: **1-800-429-5001**

9 am–6 pm PST, Mon–Fri

At all other times, please contact your healthcare professional.

### ♥ EMC Guidance

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment TMB-2080-K including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
2. Guidance and manufacturer's declaration-electromagnetic emissions and Immunity.

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class [ B ]
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Comply



Table 2

Guidance and manufacturer's declaration – electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV for power supply lines Not applicable 100 kHz repetition frequency
Surge IEC61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity								
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)	Compliance level (V/m)
	385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27	27
	450	430-470	GMRS 460, FRS 460	FM $\pm$ 5kHz deviation 1kHz sine	2	0.3	28	28
	710	704-787	LTE Band 13,17	Pulse modulation 217Hz	0.2	0.3	9	9
	745							
	780							
	810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	0.2	0.3	28	28
	870							
	930							
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS	Pulse modulation 217Hz	2	0.3	28	28
	1845							
	1970							
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9
	5500							
	5785							

## ♥ FCC Statement

### FCC ID: OU9TMB2080-B

**This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:**

**(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.**

Caution: The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

# EZFAST

Distributed by: **i-SENS USA, Inc.**  
3545 Lomita Blvd., Unit D, Torrance, CA 90505 USA  
Customer Service: **1-800-429-5001**  
9 am–6 pm PST, Mon–Fri  
At all other times, please contact your healthcare professional.