


Body Analyser Scales



Specifications

Structure

Product name	Body analyser scales					
Dimension	300x300x22.2mm (Approximately)	-		Г		Set Key
Net weight	Approximately 1.28kg (Excluding the cells)					LCD Display
Display	Digital LCD	Down Key		-)	\rightarrow	Up Key
Measurement unit	Kilogram / Pound / Stone	Electrodes -	-			- Electrodes
Measurement range	150kg ; 330lb ; 23st 8lb	-	4	$\overline{\mathbf{O}}$	\ominus	
Division	0.1kg / 0.2lb					
Accuracy	50kg: ±0.2kg; 150kg: ±0.4kg; 100kg: ±0.3kg					
Working environment	Temperature: 5°C to 40°C Relative Humidity: ≤85% RH					
	Atmospheric pressure: 86kPa to 106kPa					
Storage environment	Temperature: -20°C to 60°C Relative Humidity: ≤90% R					
	Atmospheric pressure: 50kPa to 106kPaH					
Power source	3V (One CR2032 Battery)					
Turn on method	SENSE ON technology					
Auto-off	The scale will turn off after about 10 seconds if there is no ope	ration				
Accessories	1. One CR2032 battery / 2. User Manual					
Mode of operation	Continuous Operation					
Degree of protection	Type BF applied part					

kinetik

Kinetik Medical Devices Ltd Unit 3 Perrywood Business Park, Honeycrock Lane, Surrey RH1 5DZ www.kinetikwellbeing.com

#BeWellStayWell